


**3M** Science.  
Applied to Life.™



# Cleaning and disinfection for a healthier living

## 3M™ Surface Disinfectant Cleaner

A 3-Minute, Ready-To-Use Disinfectant, Decontaminant and Cleaner for Hospital, Medical, Surgical and Institutional Use on Inanimate Environmental and Equipment Surfaces



# 3M™ Surface Disinfectant Cleaner

Fungicidal · Bactericidal · Virucidal · Tuberculocidal

A 3-Minute, Ready-To-Use Disinfectant, Decontaminant and Cleaner for Hospital, Medical, Surgical and Institutional Use on Inanimate Environmental and Equipment Surfaces



**DESCRIPTION:** 3M™ Surface Disinfectant Cleaner is a ready-to-use cleaner disinfectant that provides effective disinfection in hospitals, surgical facilities, neonatal and infant care units, child care facilities, plasmapheresis and blood centres, GI/pulmonary endoscopy labs, dialysis units and the many other areas where strict control of cross contamination is critical.

- **3M™ Surface Disinfectant Cleaner's** surfactant rich formula will quickly loosen and clean both wet and dried bodily fluids and other organic and non-organic soil from surfaces. It safely and effectively disinfects environmental and equipment surfaces such as operating room lights and tables, slit lamps, I.V. pumps, chest pumps, EMS equipment, examination room fixtures, dialysis machine surfaces, infant care equipment and incubators, oxygen hoods, respiratory therapy and anesthesia equipment surfaces, laboratory surfaces as well as other contaminated items.
- **3M™ Surface Disinfectant Cleaner** will pre-clean blood-soiled surgical instruments prior to handling during the terminal cleaning and reprocessing steps by medical device reprocessing technicians and nursing personnel. To use for this purpose, spray instruments until thoroughly drenched, allow to remain wet for 3 minutes and rinse under cool running water. Follow with an appropriate terminal cleaning, wrapping or packaging and sterilization process.

3M™ Surface Disinfectant solution provides broad spectrum accelerated quick kill of pathogenic microorganisms including:

**Gram Positive and Gram Negative Bacteria:**

- Acinetobacter baumannii*
- Mycobacterium Tuberculosis* (TB)
- Staphylococcus aureus*
- Pseudomonas aeruginosa*
- Salmonella choleraesuis* (enterica)
- Vancomycin-resistant *Enterococcus* (VRE)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)

**Viruses:**

- Both Non-enveloped and Enveloped viruses including:
- Rotavirus
  - Hepatitis B Virus (HBV)
  - Hepatitis C Virus (HCV)
  - Influenza Virus
  - Human Immunodeficiency Virus (HIV)
  - Infectious Rhinovirus
  - Herpes Simplex Virus
  - Poliovirus
  - Avian Influenza Virus (Bird Flu)

**Pathogenic Fungus:**

- Trichophyton mentagrophytes*

**3M™ Surface Disinfectant Cleaner** is safe for use on inanimate surfaces including plastics (polyvinyl chloride, polypropylene, polycarbonate and polystyrene), stainless steel, brass, aluminum, chrome, acrylic, vinyl, glass, Plexiglas®, and baked-on painted surfaces. Leaves no hazardous chemical residue when used as directed. Contains no dangerous phenols, caustic hydrogen peroxide, toxic chlorine or staining dyes. Requires no special handling, ground transportation or special storage precautions.

Unlike many other spray bottles, the specially designed 3M™ Surface Disinfectant Cleaner sprayer head produces large liquid droplets without producing fine aerosolized (airborne mist) particles that may be dispersed into the air and subsequently aspirated.

3M™ Surface Disinfectant meets and exceeds blood-borne pathogen requirements for disinfection in healthcare facilities and is formulated following strict Good Manufacturing Practices.

**APPLICATIONS:**

**3M™ Surface Disinfectant Cleaner** is a very effective:

- Surface Disinfectant and Cleaner
- Instrument Decontaminant Solution
- Ultrasonic Decontamination / Cleaning Solution

## Ordering Information

Product Description	Catalogue No.	Packaging
<b>3M™ Surface Disinfectant Cleaner (Solution) DIN 02354373</b>		
709 mL trigger spray bottle	HSD-024	12/case
3.78 L pour bottle	HSD-128	4/case
9.46 L with spigot	HSD-320	2/case
<b>3M™ Surface Disinfectant Cleaner (Wipes) DIN 02354381</b>		
6" x 10" Wipes Saturated with 3M™ Surface Disinfectant Solution packed in a Dispenser Canister	HSW-100	100 wipes/canister, 6 canisters/case
6" x 10" 3M™ Surface Disinfectant Cleaner Wipes in a peel pack	HSW-25	25 wipes/pack, 24 packs/case



Available in Canada from:

**3M Infection Prevention Division**  
**3M Canada**  
 300 Tartan Drive  
 London, Ontario N5V 4M9  
 1-800-364-3577  
 www.3M.ca

**3M Health Care**  
 2510 Conway Avenue  
 St. Paul, MN 55144-1000  
 USA  
 1-800-228-3957  
 www.3M.com

3M and 3M Science. Applied to Life. are trademarks of 3M. Used under license in Canada. Plexiglas is a registered trademark of Altuglas International. © 2017, 3M. All rights reserved. 1702-07897 iii E